
Ready to use bottle with screw-on cap, content 120ml. Clear, indelible graduation (2ml). Suitable for both single and multiple use.

Intended use:

- ☑ Collecting breast milk when used with fitting breast pump set.
- ☑ Storing of collected breast milk in refrigerator (short term) or freezer (for longer period).
- ☑ Preparing infant formula milk.
- ☐ Feeding infants when used with suitable teats or accessories. (We advise to use Materni teats & accessories)

Technical features:

- Both bottle and cap are made of food approved polypropylene. (No BPA or phthalates)
- Clear and indelible graduation for easy measuring of the pumped or fed milk volume.
- For hygienic reasons, all botlles are sealed with a srew-on cap.
- Clean by manufactoring (CBM), ready to use, no prior cleaning required.
- Suitable for bottle warmer and microwave.
- Suitable for single use (disposable) and reuse*.
- Dishwasher and freezer safe.
- Suitable for autoclave.



Available Products: Bottle 120ml with cap.

Product ID (Articlcode)	Screw Version		CBM	Qty / Box
	EU	FR		(# bottles)
MEPA0108	Χ		Χ	35
MEPA0208		Χ	Χ	35

^{* 1} box containing 6 packs of 35 pieces / total 4.200 bottles per palet.

Regulatory features:

- ☐ All used materials are food approved according regulations (EU) n° 1935:2004 and (EU) n°10-2011.
- ☑ Products are compliant to standard EN14350:2004 (1&2).
- ☐ Fully automated production (incl. packaging) under laminar flow. Cleanroom ISO Class 6 acc. to ISO14644 standard.
- ☑ Monitoring of Bioburden and microbiological testing for absence pathogene bacteria cfr. reg. (EU) n°2073:2005.

* Remark - Warning:

Single use - Multiple use.

Materni (CBM) bottles are suitable for both single and multiple use.

In case of multiple use, after first use, proper cleaning, desinfection or sterilization of the product is required to reduce the risk of (cross)contamination.

Materni strongly recommends to avoid the bottles being used by multiple infants.

For more information, instructions and safeguards, please read the instruction for use notice.

In case of multiple use, Materni declines all responsability regarding the microbiological conditions of the product.



Breastfeeding is best for babies. Before changing to formula milks, it's recommended to always consult a professional for advise.

Materni bottles, teats and accessories are suitable for professional & personal use.



